



**Sigma 8 KS Zentrifuge für TPP - TubeSpin® Bioreaktor 600/
 Sigma 8 KS Centrifuge for Use TPP TubeSpin® Bioreactor 600**

(in andere Sigmageräte passt der Bioreaktor 600 nicht hinein! / *only Sigma fit for Bioreactor 600!*)



Bild / Picture	Detail	Sigma Order #
	Rotor	• 11806
	Kapazität / Capacity	• 12 x Bioreactor 600
	Becher / Cup	• 6 x 90858
	Adapter	• 6 x 90860
	Rotor	• 11805
	Kapazität / Capacity	• 6 x Bioreactor 600
	Becher / Cup	• 1s = 6 x 13850
	Adapter	• 3s = 6 x 13851 + 3s = 6 x 17682

Abmessungen	Measurements	87600
Länge mm	Length mm	181
Durchmesser mm	Diameter mm	98
(g) max	(gx) max	3'500
Form	Form	Konisch/ conical