


# EX Pell and EX PellPLUS

## COMPATIBILITY CHART

For tips' compatibility with other pipette brands, please contact [info@capp.dk](mailto:info@capp.dk)



5130010C 5030010C 5030012C 5030015C 5030020C 5030035C	5030029C 5030030C 5130030C 5131030C	5130040C 5030040C 5030042C 5030045C 5030050C 5030055C	5030059C 5030060C 5130060C 5131060C	5130041C 5030041C 5030046C 5030051C	5030061C 5130061C 5131061C 5030063C	5030062C 5130062C 5131062C	5030064C 5030066C 5130066C 5131066C
--	--	--	--	--	--	----------------------------------	--


CAPP Bravo, ecopipette and CAPP Aero single channel pipette	0.2-2µL 0.5-10µL	0.2-2µL 0.5-10µL	0.2-2µL 0.5-10µL	0.2-2µL 0.5-10µL	0.2-2µL 0.5-10µL	0.2-2µL 0.5-10µL	2-20µL	5-50µL 10-100µL
CAPP electronic pipettes and CAPP Aero multichannel	0.5-10µL 2-20µL	0.5-10µL 2-20µL	0.5-10µL 2-20µL	0.5-10µL 2-20µL	0.5-10µL 2-20µL	0.5-10µL 2-20µL	x	5-50µL 10-100µL
Biohit eline, mline, Proline Plus	0.1-3µL 0.5-10µL	0.1-3µL 0.5-10µL	0.1-3µL 0.5-10µL	0.1-3µL 0.5-10µL	0.1-3µL 0.5-10µL	0.1-3µL 0.5-10µL	2-20µL	2-20µL 10-100µL
Brand Transferpette, Transferpette S, Transferpette Electronic	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	2-20µL	x
Eppendorf Reference, Research Plus	0.1-2.5µL 0.5-10µL	0.5-10µL	0.1-2.5µL 0.5-10µL	0.1-2.5µL 0.5-10µL	0.1-2.5µL 0.5-10µL	0.1-2.5µL 0.5-10µL	2-20µL	2-20µL 10-100µL
Finnpipette Digital	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	0.5-10µL	2-20µL	5-50µL
Gilson Pipetman	0.2-2µL 1-10µL	0.2-2µL 1-10µL	0.2-2µL 1-10µL	0.2-2µL 1-10µL	0.2-2µL 1-10µL	0.2-2µL 1-10µL	2-20µL	2-20µL 10-100µL
Rainin Classic, universal shaft (non LTS)	0.1-2µL 0.5-10µL	0.1-2µL 0.5-10µL	0.2-2µL 1-10µL	0.1-2µL 0.5-10µL	0.2-2µL 1-10µL	0.1-2µL 0.5-10µL	2-20µL	2-20µL 10-100µL
VWR, Labnet, HTL mechanical pipettes	0.1-2µL 0.5-10µL	0.1-2µL 0.5-10µL	0.2-2µL 1-10µL	0.1-2µL 0.5-10µL	0.2-2µL 1-10µL	0.1-2µL 0.5-10µL	2-20µL	2-20µL 10-100µL

© 2015 CAPP ApS. All rights reserved. CAPP Bravo, ecopipette and CAPP Aero are trademarks of CAPP ApS. eline, mline and Proline Plus are trademarks of Sartorius Biohit Liquid Handling Oy. Transferpette is a trademark of BrandTech Scientific Inc. Eppendorf Reference and Eppendorf Research Plus are registered trademarks of Eppendorf AG. Finn pipette is a trademark of Thermo Fisher Scientific Inc. Pipetman is a trademark of Gilson, Inc. Rainin Classic is a trademark of Mettler-Toledo International, Inc. WJR brand is a trademark of WJR International, Inc. Labmate is a trademark of HTLPZ HTL SA. Biopette is a trademark of Labnet International, Inc.

# EX Pell and EX Pell PLUS

## COMPATIBILITY CHART

For tips' compatibility with other pipette brands, please contact [info@capp.dk](mailto:info@capp.dk)



5130070C 5030070C 5030072C 4130075C 5030075C 5030080C 5030095C	5030089C 5030090C 5130090C 5131090C	5130100C 5130105C 5130110C	5030114C 5030115C	5130130C 5130132C 4130135C 5130135C 5130140C 5130155C	5130149C 5130150C 5131150C	5130120C 5130121C 5130122C	5130119C 5130123C	5130160C 5130165C	5130190C 5130195C
2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL
5-50µL 10-100µL 20-200µL	5-50µL 10-100µL 20-200µL	5-50µL 10-100µL 20-200µL 30-300µL	5-50µL 10-100µL 20-200µL 30-300µL	100-1200µL	100-1200µL	100-1200µL	100-1200µL	500-5000µL	x
2-20µL 10-100µL 20-200µL	10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	1-10mL
2-20µL 10-100µL 20-200µL	20-200µL	20-200µL 15-300µL	20-200µL 15-300µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	x
2-20µL 10-100µL 20-200µL	20-200µL	10-100µL 20-200µL 30-300µL	10-100µL 20-200µL 30-300µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	1-10mL
2-20µL 5-50µL 20-200µL	5-50µL 20-200µL	5-50µL 20-200µL 30-300µL	5-50µL 20-200µL 30-300µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	x	x
2-20µL 20-100µL 50-200µL	2-20µL 10-100µL 20-200µL	2-20µL 20-100µL 50-200µL	2-20µL 20-100µL 50-200µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	1-5mL	1-10mL
2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL
2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL

© 2015 CAPP ApS. All rights reserved, CAPP Bravo, ecopipette and CAPP Aero are trademarks of CAPP ApS. eline, mline and Proline Plus are trademarks of Sartorius Biohit Liquid Handling Oy. Transferpette is a trademark of BrandTech Scientific Inc. Eppendorf Reference and Eppendorf Research Plus are registered trademarks of Eppendorf AG. Finn pipette is a trademark of Thermo Fisher Scientific Inc. Pipetman is a trademark of Gilson, Inc. Rainin Classic is a trademark of Mettler-Toledo International, Inc. WJR brand is a trademark of WJR International, Inc. Labmate is a trademark of HTLPZ HTL SA. Biopette is a trademark of Labnet International, Inc.



