



Expell 384 Plates Compatibility Chart
Thermocycler Compatibility for 5100513C
ref. PIP_DS_5100513C_22082018

| | | | | |
|----------------------|---------------------------------------|-------------------------------|---|------------|
| Agilent | SureCycler 8800 thermal cycler | 384 microplate | ✓ | Compatible |
| Analytik-Jena | FlexCycler | Monoblock 384 | ✓ | Compatible |
| Analytik-Jena | SpeedCycler 96 | 384-well | ✓ | Compatible |
| Apollo | ATC401 | 384 microplate | ✓ | Compatible |
| Apollo | ATC701 | 384 microplate | ✓ | Compatible |
| Apollo | ATC801 | 384 microplate | ✓ | Compatible |
| Applied Biosystems | 3500xL Genetic Analyzer | 24-capillary | ✓ | Compatible |
| Applied Biosystems | 7900HT Fast Real-Time PCR System | 384-well | ✓ | Compatible |
| Applied Biosystems | GeneAmp PCR System 9700 | Dual 384-Well 384 microplate | ✓ | Compatible |
| Applied Biosystems | QuantStudio™ 12k Flex | 384-well | ✓ | Compatible |
| Applied Biosystems | QuantStudio™ 5 Real-Time PCR System | 384-well | ✓ | Compatible |
| Applied Biosystems | QuantStudio™ 6 Flex | 384-well | ✓ | Compatible |
| Applied Biosystems | QuantStudio™ 7 Flex | 384-well | ✓ | Compatible |
| Applied Biosystems | Veriti™ 384-Well Thermal Cycler | 384 microplate | ✓ | Compatible |
| Applied Biosystems | Veriti™ Dx 384-Well Thermal Cycler | 384 microplate | ✓ | Compatible |
| Applied Biosystems | ViiA™ 7 Real-Time PCR System | 384-well | ✓ | Compatible |
| Bioer | GenePro B-384G | 384 microplate | ✓ | Compatible |
| Bioer | XPcycler TC-XP-F | 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | C1000 | 384-Well | ✓ | Compatible |
| Bio-Rad Laboratories | C1000 Touch | 384-Well | ✓ | Compatible |
| Bio-Rad Laboratories | CFX384 Real-Time PCR Detection System | 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | DNA Engine | 384 Alpha unit 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | DNA Engine Dyad | 384 Alpha unit 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | DNA Engine Dyad Disciple | 384 Alpha unit 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | DNA Engine Tetrad 2 | 384 Alpha unit 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | iCycler | 384 microplate | ✓ | Compatible |
| Bio-Rad Laboratories | S1000 | 384-Well | ✓ | Compatible |

| | | | | |
|--|---------------------------------|-------------------------------|---|------------|
| Boeco | TC-PRO 384G | 384 microplate | ✓ | Compatible |
| Bulldog Bio | GenePro™ Thermal Cycler | 384-well | ✓ | Compatible |
| Bulldog Bio | GeneTouch™ Thermal Cycler | 384-well | ✓ | Compatible |
| Corbett Research | Palm-Cycler | 384 microplate | ✓ | Compatible |
| eENZYME | GenePro Thermal Cycler | 384 microplate | ✓ | Compatible |
| Eppendorf | Mastercycler ep | 384 microplate | ✓ | Compatible |
| Eppendorf | Mastercycler pro 384 | 384 microplate | ✓ | Compatible |
| Eppendorf | Thermomixer C | 384-well | ✓ | Compatible |
| Esco | Swift™ Maxi | 384 microplate | ✓ | Compatible |
| Esco | Swift™ MaxPro | 384-wells | ✓ | Compatible |
| EuroClone | peqSTAR 384 HPL | 384 microplate | ✓ | Compatible |
| EuroClone | peqSTAR 384 ML | 384 microplate | ✓ | Compatible |
| Finnzymes Instruments | TCA0384 Arktik | 384 | ✓ | Compatible |
| Gene Technologies Ltd. | G-Storm GS1 | 384 microplate | ✓ | Compatible |
| Gene Technologies Ltd. | G-Storm GS2 | 384 microplate | ✓ | Compatible |
| Gene Technologies Ltd. | G-Storm GS4 | 384 microplate | ✓ | Compatible |
| Gene Technologies Ltd. | G-Storm GSX | 384 microplate | ✓ | Compatible |
| Hangzhou Jingle Scientific Instrument Co.,LTD. | KF960-D | 384-well | ✓ | Compatible |
| LongGene Scientific Instruments Co.,Ltd. | MyGene™ PCR | 384-wells | ✓ | Compatible |
| MJ Research Inc. | PTC-200 DNA Engine | 384 Alpha unit 384 microplate | ✓ | Compatible |
| MJ Research Inc. | PTC-220 DNA Engine Dyad | 384 Alpha unit 384 microplate | ✓ | Compatible |
| MJ Research Inc. | PTC-221 DNA Engine Disciple | 384 Alpha unit 384 microplate | ✓ | Compatible |
| MJ Research Inc. | PTC-225 DNA Engine Tetrad | 384 Alpha unit 384 microplate | ✓ | Compatible |
| MJ Research Inc. | PTC-240 DNA Engine Tetrad 2 | 384 Alpha unit 384 microplate | ✓ | Compatible |
| MWG-Biotech | Primus Dualblock PCR-System | 384 microplate | ✓ | Compatible |
| MWG-Biotech | Primus HT Multiblock PCR-System | 384 microplate | ✓ | Compatible |
| MWG-Biotech | Primus PCR-System | 384 microplate | ✓ | Compatible |
| PCRMax | AlphaCycler AC-1 | 384-well | ✓ | Compatible |
| PCRMax | AlphaCycler AC-4 | 384-well | ✓ | Compatible |
| Peqlab | peqSTAR 384 ML | 384 microplate | ✓ | Compatible |
| Peqlab | peqSTAR 384 MPL | 384 microplate | ✓ | Compatible |
| Peqlab | Primus 96 Advanced | 384 microplate | ✓ | Compatible |
| Peqlab | Primus 96 Advanced HT2X | 384 microplate | ✓ | Compatible |

| | | | | |
|-------------------------|--|-----------------|---|------------|
| Peqlab | Primus 96 Advanced HT4X | 384 microplate | ✓ | Compatible |
| Roche | LightCycler ® 480 Real-Time PCR System | 384-well | ✓ | Compatible |
| Scinics | EZ-Cycler MP-384 | 384 microplate | ✓ | Compatible |
| SensoQuest | Labcycler | Thermoblock 384 | ✓ | Compatible |
| Techne | Flexigene | 384 microplate | ✓ | Compatible |
| Techne | Genius | 384 microplate | ✓ | Compatible |
| Techne | Prime thermal cycler | 384 microplate | ✓ | Compatible |
| Techne | TC-4000 | 384 microplate | ✓ | Compatible |
| Techne | TC-412 | 384 microplate | ✓ | Compatible |
| Techne | TC-5000 | 384 microplate | ✓ | Compatible |
| Techne | TC-512 | 384 microplate | ✓ | Compatible |
| Techne | TC-Plus | 384 microplate | ✓ | Compatible |
| Thermo Hybaid | MultiBlock System | 384 microplate | ✓ | Compatible |
| Thermo Hybaid | PCR Express | 384 microplate | ✓ | Compatible |
| Thermo Hybaid | Px2 | 384 microplate | ✓ | Compatible |
| ThermoFisher Scientific | Arktik Thermal Cycler | 384-well | ✓ | Compatible |
| VWR | DuoCycler | 384 microplate | ✓ | Compatible |
| VWR | DuoCycler fast | 384 microplate | ✓ | Compatible |
| VWR | QuattroCycler | 384 microplate | ✓ | Compatible |
| VWR | QuattroCycler fast | 384 microplate | ✓ | Compatible |
| VWR | UnoCycler | 384 microplate | ✓ | Compatible |
| VWR | UnoCycler fast | 384 microplate | ✓ | Compatible |
| Whatman Biometra | T1 Thermocycler | 384 microplate | ✓ | Compatible |
| Whatman Biometra | TProfessional | 384 microplate | ✓ | Compatible |
| Whatman Biometra | TRobot | 384 microplate | ✓ | Compatible |
| Whatman Biometra | UNO II | 384 microplate | ✓ | Compatible |